

# TEST SERVICE

Accredited Calibration and Verification Service  
for Balances, Test Weights, Force Measurement and more



Deutsche  
Akkreditierungsstelle  
D-K-19408-01-00

Accredited calibrations according to  
DIN EN ISO/IEC 17025:2018 norms for balances,  
weights, force, volumes of solid bodies, densities  
of solid bodies, temperature, humidity.

**kern-lab.com**



# 2026

# CONTENTS

The heart of calibration and verification .....	3
All you need to know about calibration .....	4
KERN test services at a glance .....	5
The balance .....	6
Calibration and verification prices for electronic balances .....	11
Equipment qualification .....	12
The test weight .....	13
The force gauge .....	18
Factory calibration/temperature and relative humidity .....	20
Digital calibration certificate (DCC) .....	22
KERN glossary (contains the key word which are marked with ►) .....	23

## Your contact persons at the KERN accredited calibration laboratory

**Test services sales** +49 7433 9933-196 // [recalibration@kern-sohn.com](mailto:recalibration@kern-sohn.com)



Team leader  
**Lars Wagner**  
+49 7433 9933-255  
[lars.wagner@kern-sohn.com](mailto:lars.wagner@kern-sohn.com)



Deputy team leader  
**Waldemar Fleitling**  
+49 7433 9933-163  
[fleitling@kern-sohn.com](mailto:fleitling@kern-sohn.com)



Sales Manager  
**Tülin Lök**  
+49 7433 9933-148  
[tuelin.loek@kern-sohn.com](mailto:tuelin.loek@kern-sohn.com)



Backoffice  
**Katharina Heise**  
+49 7433 9933-172  
[katharina.heise@kern-sohn.com](mailto:katharina.heise@kern-sohn.com)



Backoffice  
**Annalena Wuhler**  
+49 7433 9933-217  
[annalena.wuhler@kern-sohn.com](mailto:annalena.wuhler@kern-sohn.com)



Backoffice  
**Balbina Pietrzak**  
+49 7433 9933-210  
[balbina.pietrzak@kern-sohn.com](mailto:balbina.pietrzak@kern-sohn.com)



Backoffice  
**Alexander Jonitz**  
+49 7433 9933-183  
[alexander.jonitz@kern-sohn.com](mailto:alexander.jonitz@kern-sohn.com)

**Test services on-site** +49 7433 9933-400 // [testservices-onsite@kern-sohn.com](mailto:testservices-onsite@kern-sohn.com)



Team leader  
**Karl-Richard Fuchs**  
+49 7433 9933-136  
[fuchs@kern-sohn.com](mailto:fuchs@kern-sohn.com)



Deputy team leader  
**Lars Wagner**  
+49 7433 9933-255  
[lars.wagner@kern-sohn.com](mailto:lars.wagner@kern-sohn.com)



Backoffice  
**Tanja Michailidis**  
+49 7433 9933-134  
[tanja.michailidis@kern-sohn.com](mailto:tanja.michailidis@kern-sohn.com)

**Calibration partnership** +49 7433 9933-345 // [calpartner@kern-sohn.com](mailto:calpartner@kern-sohn.com)



Team leader  
**Michael Stark**  
+49 7433 9933-164  
[stark@kern-sohn.com](mailto:stark@kern-sohn.com)



Backoffice  
**Martin Goltz**  
+49 7433 9933-151  
[goltz@kern-sohn.com](mailto:goltz@kern-sohn.com)

# The heart of calibration and verification\*<sup>1</sup>

## The principle

Every electronic measuring device will only give correct results if it is checked regularly, i.e. calibrated correctly and adjusted when required. An electronic balance, test weight or another measuring device is only a reliable measuring and checking tool if it is calibrated and this calibration is documented as part of a quality procedure.

Calibration accredited by DAkkS (DAkkS = German calibration service) documents traceability to the national standard and this then meets the standard requirements of QM systems. Accredited calibrations are valid internationally.

## Calibration of measuring devices


Measuring “correctly” is of elementary significance, as it is not unusual for inaccurate or “wrong” measurements to have expensive economic consequences. Calibration or establishing the accuracy of checking equipment is carried out by accredited laboratories throughout the world in accordance with the DIN EN ISO/IEC 17025 standard. On an international level, the EA (European co-operation for Accreditation) and ILAC (International Laboratory Accreditation Cooperation) monitor the upholding of the highest quality standards. In Germany this is carried out by DAkkS (German accreditation point).

## What does calibration mean?

Determining and documenting the deviation from true, actual measure value of the value displayed by a measuring device or of the value given by checking equipment.



Recalibration information (optional)



DAkkS calibration mark

## When is an accredited calibration required?

Accredited calibration is always necessary, when checking equipment is to be used in a QM process (e.g. in accordance with ISO 9000ff, TS 16949, VDA, FDA, GLP, GMP, ...). The operator controls the use of checking equipment and periodic recalibration time intervals themselves. Calibration certificates with accreditation symbol are recognised internationally.

## Deutsche Akkreditierungsstelle (DAkkS)

The German accreditation body (DAkkS) is the successor to the German calibration service (DKD) in terms of accreditation systems. On the basis of EC regulation no. 765/2008, the accreditation point of the German calibration service (DKD) was transferred to the German accreditation point (DAkkS) with effect from 17.12.2009.

## Who needs a calibration certificate with accreditation symbol?

In the context of standard requirements for monitoring checking equipment, every company with a Quality Management system is obliged to test and document its measuring equipment at regular intervals. A calibration certificate with accreditation symbol fulfils this obligation.

## The KERN calibration laboratory (D-K-19408-01-00)

KERN has a highly-automated calibration laboratory with accreditation to DIN EN ISO/IEC 17025 in the field of balances, test weights, force measurement, humidity and temperature. By using the most modern calibration technology with high-end calibration robots in fully air-conditioned laboratories, the measurement uncertainty and process times are reduced to a minimum, and also the quality of the calibration is increased. As an accredited and certified calibration service provider with decades of experience, we offer you an extensive range of services, which will leave no demand unfulfilled. The accreditation applies to the extent specified in the appendix to the certificate D-K-19408-01-00.

## Calibration or verification

**Accredited calibration** is possible for every balance in perfect condition. Accredited calibration is a private service for ensuring high quality requirements according to DIN EN ISO 9000ff and other standards, e.g. in production and research. Verifying is only possible for balances with certificate of compliance marked with a green **M**.

More interesting facts at: [www.kern-lab.com](http://www.kern-lab.com)

\*<sup>1</sup> The “initial verification” for new balances is called conformity assessment according to NAWID: 2014/31/EU, a verification corresponds to the “reverification”.



# All you need to know about calibration and verification\*1



## Accredited calibration

(area not regulated by law)

### Why?

Accredited calibration is always necessary when checking equipment (balance or test weight) is to be used in a QM process (e.g. to ISO 9000ff, GS 9000, TS 16949, VDA 6.1, FDA, GLP, GMP, ...)

### What?

An accredited calibration can be performed on any checking equipment.

### How?

Determination of accuracy throughout the world by a laboratory which is accredited to DIN EN ISO/IEC 17025. Traceability to internationally recognised standards. The calibration certificate with accreditation symbol confirms both the metrological characteristics of the checking equipment and the general requirements for the ► *monitoring of checking equipment* (eg. ISO 9000ff).

### Where?

Internationally recognised. This is monitored by EA (European co-operation for Accreditation) and ILAC (International Laboratory Accreditation Cooperation), and in Germany, for example the DAkkS (Deutsche AkkreditierungsStelle GmbH) – German accreditation point.

### When?

The operator controls the use of checking equipment and periodic recalibration time intervals themselves.

## Verification\*1

(area regulated by law)

### Why?

Applications with mandatory verification of balances and test weights include commercial trade when the price of a commodity is determined by weighing, the manufacture of pharmaceuticals in pharmacies, the production of pre-packaged goods in medical applications.

### What?

You can only verify balances which have a certificate of compliance and test weights which conform to ► *OIML* standards.

### How?

Testing to verification permissible error limits (for details on tolerances see page 14) to protect the consumer. When introducing balances and weights onto the market, they are subject to EU directives. The subsequent monitoring of the market is regulated at a national level, in Germany through the MessEG (Weights and Measures Act) and MessEV (Verification ordinance).

### Where?

EC Declaration of Conformity with CE marking is valid as “Initial verification” throughout Europe. Reverification and national declarations of conformity are only recognised on a national level.

### When?

The legislative body governs the use of balances and test weights as well as time intervals for re-verification. National specifications apply here.

\*1 The “initial verification” for new balances is called conformity assessment according to NAWID: 2014/31/EU, a verification corresponds to the “reverification”.



# KERN test services at a glance

## Calibration of balances inhouse (at KERN)

With the shortest calibration time in the KERN calibration laboratory of maximum 4 working days after receipt of order, this gives you almost uninterrupted use of your balances within your production process.s.

## Calibration of balances on site (at the customer)

As an option, you can have your balances calibrated on your premises. This on-site testing service is metrologically recommended, as the balance is in its field of use and can be calibrated without any possible transportation problems. Minimized downtime and personal contact with our expert are the major benefits of this service. We would be pleased to give you more information and agree a date with you.

## Calibration of weights

Here too, with its short process times, KERN is unbeatable. The most modern calibration robots calibrate your test weights with only the slightest ► *measuring uncertainty*, according to international directives of OIML R 111 and thereby ensuring a reliable weighing result. Recommended recalibration period 1 year. On-site calibration of your weights according to OIML classes M1 – M3 (10 kg – 50 kg) can also offer you an affordable alternative. We would be pleased to come to your premises and calibrate your test weights with our mobile **MACOS calibration system**.

## Calibration of force

Through the force-measurement accreditation from KERN (in Newtons), Accredited calibration of your force-measuring devices means that we can meet the highest requirements. With test stands and measuring procedures designed for this purpose, our specialists can calibrate your checking equipment to the latest test methodology in our laboratory.

## Calibration of temperature and humidity

With the shortest calibration time in the KERN calibration laboratory of maximum 4 working days after receipt of order, this gives you almost uninterrupted use of your balances within your production process.s.

## Volume determination

When calibrating every new weight in OIML class E1 you must also establish its volume. This is necessary for the correction of air buoyancy. Accredited volume determination in our laboratories is an integral part of our high-end demands.

## Reconditioning of weights

KERN gets your weights back up to standard, **regardless of the manufacturer**. Whether it is adjustment, marking, sand blasting or lacquering. The aim here is compliance and long-term stability. Special arrangements on request.

## Magnetic characteristics

By measuring sensitivity /and magnetisation KERN gains reliable evidence regarding the magnetic characteristics of your test weights. “Magnetic” weights can distort the weighing result when you are using the balance.

## Factory calibration

The testing of measuring devices for accuracy in accordance with a recognised but not accredited process without proof of metrological traceability – this is the difference when compared with accredited calibration.

## Digital Calibration Certificate (DCC)

The “Digital Calibration Certificate (DCC)” created by the PTB can be downloaded on [www.kern-lab.com/dcc](http://www.kern-lab.com/dcc). You will find further information on page 22.

The paperless calibration certificate as PDF can be downloaded at the calibration download on [www.kern-lab.com](http://www.kern-lab.com).

## Reverification service for balances and test weights

The reverification of balances and weights is regulated on a national level and can therefore only be offered for balances and weights used in Germany. For reverification kindly contact the agencies in your country.

## Database supported management of checking equipment

Information on your checking equipment which has been calibrated by us is stored in our database. In this way it is possible to make trend calculations. You therefore get an overview of long-term stability and trend behaviour of your checking equipment.

## Reminder service


The continuous recalibration of your checking equipment is an integral part of the reliable management of checking equipment. You can depend on KERN to support you, and KERN will remind you in good time when the next recalibration is due.

**This service is free of charge to you!**

## Collection and delivery service

Why not let us transport your checking equipment correctly. We will collect your checking equipment from you and then deliver it quickly and safely.

# The balance



## a) KERN on-site calibration (we visit you)

In Germany, KERN has a close-knit network of accredited calibration laboratory employees, who can carry out on-site calibration of balances up to 50 tonnes.

This on-site testing service is metrologically recommended, as your balance is in its field of use and can be calibrated without any possible transportation problems.

Lower downtime and personal contact with our expert are the major benefits of this service.

**This KERN calibration service is also independent of the brand.** Preparatory maintenance work by agreement. Prices for on-site calibration on request.

**Tell us your desired date with indication of the scales to be tested or enter your inquiry directly in our offer generator on [www.kern-lab.com](http://www.kern-lab.com). One of our calibration staff will then contact you immediately and discuss the calibration procedure with you at your premises – uncomplicated and competent.**

**i** phone: +49 7433 9933-400  
e-Mail: [testservices-onsite@kern-sohn.com](mailto:testservices-onsite@kern-sohn.com)

## THE ADVANTAGES OF USING THE KERN ON-SITE CALIBRATION:



- + Calibration on-site at your premises in the field of use
- + No risk of damage during transportation
- + Low downtime
- + Cross-brand servicing, basic inspection and adjustment by a specialist
- + You tell us when you would like us to come
- + Device training for qualified users



## Calibration of balances

Any balance will only give correct results if it is checked regularly, i.e. calibrated correctly and adjusted when required. A balance is only a reliable measuring and checking tool if it is calibrated and this calibration is documented. The issued accredited calibration certificates are proof of the metrological traceability to national and international standards, as required by the DIN EN ISO 9000 and DIN EN ISO/IEC 17025 standards, amongst others. KERN recommends a recalibration period of one year. The standard does not give a defined recalibration period. KERN recommends that, with intensive (daily) use, you recalibrate your balance every 6 months and with normal (weekly) use, every 12 months.



### THE ADVANTAGES OF USING THE KERN IN-HOUSE CALIBRATION:



- + Short calibration time: Test time in the laboratory is only four working days
- + Competence: Calibration laboratory, which complies with the highest standards in the area of metrology
- + Management of the recalibration calendar for your individual measuring instrument is possible
- + Cross-brand service: Measuring devices from any manufacturer can be calibrated independently
- + Repair: Any necessary repairs can be carried out immediately, if you wish



### b) Calibration at the KERN factory (you send your balance to us)

Recommended for new devices and for balances which can be affordably transported, as then there is no need for us to travel to carry out the calibration on-site. Repairs can be carried out at the same time, quickly and in full.

#### The process would be as follows:

- Day 1: Send your balance to the KERN calibration laboratory in Balingen.
- Day 2 to 3: Evaluation and calibration of your balance by our specialists.
- Day 4: After positive validation, your balance is returned.

#### Recalibration

##### • Typical industrial recalibration times

may be recommended as follows:

- daily use of the measuring instrument (once or several times): Recalibration period of 6 months
- weekly use of the measuring instrument (or less often): Recalibration period of 12 months

- **Recalibration prices:** The prices for (re)calibration of balances can be found on page 11. Costs for cleaning, function testing and any necessary adjustment or for the production of special holders to carry out the calibration will be calculated separately.



## Calibration certificate with accreditation symbol for balances

### 1. "Official" document

The calibration laboratory KERN (D-K-19408-01-00) is accredited through the accreditation point of the Deutsche Akkreditierungsstelle GmbH. The calibration certificate with accreditation symbol is recognised internationally and is available in several languages.

### 2. Item to be calibrated

The calibration item as well as the type or model with serial number is documented. This means that there is no confusion and guarantees the assignment of the calibration certificate to a specific balance.

### 3. Traceability

The reference standards of the accredited laboratory are monitored in strictly defined cycles and periodically brought into line with national and thereby international standards. This is carefully documented and given on the calibration certificate. In this way the basic fundamental traceability to the national standard is ensured.

### 4. Applicant

On the very first page of the calibration certificate you will clearly see the applicant or owner of the calibrated checking equipment.

### 5. Metrological part

As well as other tests, three metrological tests are carried out during accredited calibration. These are to test repeatability, accuracy and eccentric loading. This defines the features of the balance.

### 6. Measurement uncertainty of a balance

This is determined individually for each balance according to a precisely given test method and is documented in the calibration certificate. It depends on various factors, both internal and external to the balance.

### 7. Usage accuracy

Usage accuracy gives the uncertainty when the operator uses the measuring equipment on site. This value, which is established by a mathematical equation is influenced by changes in temperature, type of use and other factors.

### 8. Minimum weight of sample

(optional; see page 9 Art.No. 969-103)

The smaller the sample weight, the larger the relative measuring uncertainty. For those responsible for weighing processes, it is important to determine the deviations which occur when establishing values of the smallest loads. Determining the minimum sample weight declares in this way, clearly the various requirements on the weighing accuracy in relation to the sample weight.

**KERN® KERN & SOHN GmbH**  
Akreditiertes Kalibrierlabor seit 1994.  
Accredited calibration laboratory since 1994.  
Ihr Partner für Kalibrierdienstleistungen, Prüfmittelmanagement und Beratung.  
Your partner for calibration services, test equipment management and support.

Mitglied im / member of the  
**Deutschen Kalibrierdienst** **DKD** **ILAC-MRA** **DAKKS**  
Deutsche Akkreditierungsstelle  
D-K-19408-01-00

Kalibriert nach / calibrated according to  
Sample 2026-01/1

Kalibriert nach / calibrated according to  
Sample 2026-01

Gegenstand / Object  
Analysewaage  
Analytical Balance

Hersteller / Manufacturer  
KERN & SOHN GmbH  
Ziegels 1  
72336 Balngen-Fronnen

Typ / Type  
ABT 120-SDM

Fabrikat/Serien-Nr.  
Serial number  
WX12345678

Auftraggeber / Customer  
Mustermann GmbH  
Musterweg 42  
12345 Musterstadt  
Deutschland

Dieser Kalibrierschein dokumentiert die metrologische Rückführbarkeit auf nationale Normale zur Darstellung der Einheiten in Übereinstimmung mit dem internationalen Einheitensystem (SI). Die DAKKS ist Unterzeichner der multilateralen Übereinkommen der Europäischen Kooperations für Akkreditierung (EA) und der International Laboratory Accreditation Cooperation (ILAC) zur gegenseitigen Anerkennung der Kalibrierscheine. Für die Einhaltung einer angemessenen Frist zur Wiederholung der Kalibrierung ist der Benutzer verantwortlich.

This calibration certificate documents the metrological traceability to national standards, which realize the units of measurement according to the International System of Units (SI). The DAKKS is a signatory to the multilateral agreements of the European Cooperation for Accreditation (EA) and the International Laboratory Accreditation Cooperation (ILAC) for the mutual recognition of calibration certificates. For the observance of a reasonable period for the repetition of the calibration, the user is responsible.

Messergebnisse: / Measurement results: 2026

Zustand #1: Ursprungszustand / as found  
State -

Temperatur: zu Beginn 22,0 °C  
Temperature at the beginning

Bemerkungen / Remarks:  
Der Kennwert der Waage wurde vor der Kalibrierung mit dem internen Justiergewicht justiert.  
Before calibration, the span was adjusted with the internal calibration weight.

1. Wiederholbarkeit / Repeatability

Messung / Measurement	Prüflast / Load	Waagenanzeige / Indication
No. 1	100 g	100.0002 g
No. 2	100 g	100.0003 g
No. 3	100 g	100.0004 g
No. 4	100 g	100.0004 g
No. 5	100 g	100.0004 g

Standardabweichung:  $s = 0.00009$  g  
Standard deviation

2. Außermittige Belastung / Eccentricity

Position / Position	Prüflast / Load	Waagenanzeige / Indication
No. 1	50 g	50.0001 g
No. 2	50 g	50.0001 g
No. 3	50 g	50.0001 g
No. 4	50 g	50.0001 g
No. 5	50 g	50.0002 g

Messunsicherheit / Measurement uncertainty 2026

Angaben ist die erweiterte Messunsicherheit, die sich aus der Standardmessunsicherheit durch Multiplikation mit dem Erweiterungsfaktor  $k$  ergibt. Sie wurde gemäß EA-402:2022 und EURAMET-19v4.0 ermittelt. Der Wert der Messgröße liegt im Regelfall mit einer Wahrscheinlichkeit von mindestens 95% im zugeordneten Wertebereich. Die Ergebnisse gelten nur für den kalibrierten Gegenstand in Zustand und unter den Bedingungen zum Zeitpunkt der Kalibrierung. Ein Anteil für die Langzeitstabilität des Kalibriergegenstandes ist nicht enthalten.

The expanded measuring uncertainty is calculated by multiplication of the standard measuring uncertainty with the coverage factor  $k$ . It was determined according to EA-402:2022 and EURAMET-19v4.0. The value of the test result is normally with a probability of at least 95% within the assigned value interval. The results apply only to the calibrated item in the condition and under the conditions at the time of calibration. A proportion for the long-term stability of the calibration item is not included.

Zustand / State #1 (Ursprungszustand) / as found, - 1 -

Prüflast / Load	Abweichung / Error	Erweiterungs-faktor $k$ / Coverage factor	Unsicherheit / Uncertainty	relative Unsicherheit / Rel. uncertainty
20 g	0.0001 g	2.27	0.00026 g	0.00129 %
50 g	0.0002 g	2.18	0.00028 g	0.00054 %
70 g	0.0003 g	2.05	0.00035 g	0.00049 %
100 g	0.0004 g	2.06	0.00034 g	0.00033 %
120 g	0.0005 g	2.02	0.00043 g	0.00035 %

Darstellung im Diagramm / Representation as chart

Verwendungsgenauigkeit / Total usage accuracy

Bei der Verwendung der Waage erhöht sich die Messunsicherheit aufgrund verschiedener Einflüsse. Unter Annahme der gleichen Umgebungsbedingungen (z. B. Windzug, Erschütterungen) wie bei der Kalibrierung und geschätzten Raumtemperaturschwankungen von 1 K bei einem geschätzten Temperaturkoeffizienten von 1 ppm/K ergibt sich die unten genannte Verwendungsgenauigkeit gemäß EURAMET-19v4.0. Dabei sind Abweichungen der Waage bewertungslos, die Anzeige der Waage muss also nicht korrigiert werden. Es wird davon ausgegangen, dass die Waage regelmäßig justiert wird.

Several effects increase the measuring uncertainty of the balance at ultramicro. Assuming the same environmental conditions as at calibration time with an estimated room temperature variance of 1 K and an estimated temperature-coefficient of 1 ppm/K, the following usage accuracy is calculated according to EURAMET-19v4.0. The determined errors of indication were considered, so no correction needs to be applied. It is assumed that the balance will be adjusted regularly.

$G = 0.00009 \text{ g} + 1.03 \cdot 10^{-4} \cdot m_e$

$m_e$  = Nettoanzeige bei zunehmender Belastung  
net display with increasing load

Diagramm der Verwendungsgenauigkeit / Graph of usage accuracy:

— Prozessgenauigkeit / Process accuracy  
— Messunsicherheit / Uncertainty

rel. Meßunsicherheit

Mindesteinwaage / Minimum weight of sample

Nettobelastung

Calibration certificate with accreditation symbol for balances (extract)

## Minimum weight of sample (in use)

**What is the lightest item you can weigh on your balance, while still achieving accurate and reliable weighing results? What exactly is the limit?**

The KERN minimum sample weight protocol accounts for the established minimum sample weight of your balance and its location of installation and use with the relative ► **measuring uncertainty**. With various safety coefficients and required weighing accuracy (process accuracy), depending on standard or quality-related requirements on the balance being used.

The higher the selected safety coefficient, the higher the safety when using the balance in a particular process. Typical perturbations when using the balance e.g. small fluctuations in temperature are taken into account. In easily predictable conditions in a professional environment of use, KERN recommends a safety coefficient of 3. For critical processes, a correspondingly higher factor should be selected. The minimum sample weight protocol contains a diagram as well as a table, from which you can ascertain the minimum sample weight for your balance, depending on the process.


### Safety coefficients and required process accuracy for the minimum sample weight:

#### Usage accuracy

As the calibration of a balance is a snapshot, there must be a statement relating to how the measuring instrument behaves in operation in the location of use, or to what extent the measuring uncertainty may change. This is because with daily use, the measuring uncertainty of a balance increases due to various influences. These influences must be recorded and rated ...

#### ... and this is how:

By accepting that the same environmental conditions (e.g. draughts, vibrations, ...) as they were at the point of calibration are present at the location of the balance and estimated room temperature fluctuations of X Kelvin (°C) with a temperature coefficient assigned in the balance (in ppm/K), the result is a particular accuracy of use. You can establish this accuracy of use in accordance with EURAMET/cg-18.



Example:

Balance with 220 g.

At 82.5 g the usage accuracy

is 0.0005 g.

$\hat{=}$  0.000606 %

Required process accuracy	Safety coefficient			
	1	3	5	10
0,1 %	0,0985 g	0,2983 g	0,5021 g	1,0297 g
0,2 %	0,0491 g	0,1480 g	0,2480 g	0,5021 g
0,5 %	0,0196 g	0,0590 g	0,0985 g	0,1979 g
<b>1,0 %</b>	<b>0,0098 g</b>	<b>0,0294 g</b>	<b>0,0491 g</b>	<b>0,0985 g</b>
2,0 %	0,0049 g	0,0147 g	0,0245 g	0,0491 g
5,0 %	0,0020 g	0,0059 g	0,0098 g	0,0196 g
10,0 %	0,0010 g	0,0029 g	0,0049 g	0,0098 g

# Adjustment at the location of installation

## Why?

Adjustment at the location of installation is necessary, as the measuring results of balances depend on the local gravitational force (gravitational acceleration) and therefore depend on the location of use. KERN can carry this out just before shipping at the factor, individually to suit the location of installation.

## What are the advantages of carrying out adjustment at the location of installation?

- The balance gives reliable measurement results at the location of installation.
- No time-consuming on-site adjustment necessary.
- You do not need a Service Engineer or any additional weights.
- The balance is ready for immediate use.

For adjustment to the location of installation you need the value for gravitational acceleration at the location of installation, which KERN can calculate using the point of use. The procedure is suitable for balances with a resolution of <60,000 d. For higher resolutions we recommend a balance with an internal adjusting weight or adjustment with a calibrated adjusting weight at the location of installation.

## Pricing table for adjustment at the location of installation

Area	KERN	Price
[Max] ≤ 5 kg	961-247	45,-
[Max] > 5 – 50 kg	961-248	55,-
[Max] > 50 – 350 kg	961-249	65,-
[Max] > 350 – 1500 kg	961-250	104,-
[Max] > 1500 – 2900 kg	961-251	138,-
[Max] > 2900 – 6000 kg	961-252	275,-
[Max] > 6000 – 12000 kg	961-253	315,-

**KERN**  
CALIBRATION

**KERN & SOHN GmbH**  
Kalibrierlabor seit 1994.  
Calibration laboratory since 1994.  
Ihr Partner für Kalibrierdienstleistungen, Prüfmittelmanagement und Beratung.  
Your partner for calibration services, test equipment management and support.

**Justage auf den Aufstellungsort**  
Adjustment to the place of use

Kalibriergegenstand:  
Calibration object

Hersteller  
Manufacturer

Seriennummer  
Serial no.

Auftragsnummer  
Order No.

Auftraggeber  
Customer

Ort der Justage  
Place of adjustment

IFB 30K5DM

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

DB1234567

2026-12345678

Mustermann GmbH  
Musterstr. 1  
12345 Musterstadt  
Deutschland

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen-Frommern  
Deutschland

Die Justage auf den Aufstellungsort wurde vom Kunden erwünscht. Die Waage wurde mit rückführbaren Normalen auf die angegebene Gravitation justiert. Gegen eine weitere Justage ist die Waage nicht gesichert.

The adjustment to the place of use was requested by the customer. The balance was adjusted using weights which are traceable to the national standards. The weighing instrument is not secured against a re-adjustment.

# Certificate of conformity

**Konformitätszertifikat / Certificate of conformity**  
ausgestellt für: / issued for:  
Typ: PNJ 3000-2M Serien-Nr.: WX161234567 Inventar-Nr.: -  
Type: Serial number Inventory number

Dieses Konformitätszertifikat bescheinigt die Gültigkeit der folgenden Konformitätsaussagen auf Basis der Messergebnisse des Kalibrierscheins.

This conformity certificate documents the validity of the following statements of conformity based on the measurement results of corresponding calibration certificate:

Konformitätsaussagen: / Statements of conformity:

A) Die im Rahmen der Kalibrierung festgestellten Anzeigeabweichungen der Waage (siehe Seite 4) liegen unter den angegebenen Mess- und Umgebungsbedingungen und unter Berücksichtigung der erweiterten Messunsicherheiten (Überschreitungswahrscheinlichkeit 95%) innerhalb der Toleranz. Die angegebene Messunsicherheit berücksichtigt bereits unter anderem die Einflüsse der Wiederholbarkeit und der äußeren Belastung, weshalb eine separate Bewertung dieser Parameter nicht durchgeführt wurde.

The errors of indication determined during calibration (ref. page 4) are under given measurement and environmental conditions and considering the expanded measurement uncertainty (coverage probability 95%) within the specification. The given measurement uncertainty already includes the effects of repeatability and eccentricity. Therefore, those parameters have not been assessed separately.

Nr. No.	Prüflast Test load	Anzeige Indication	Abweichung Error	erw. Unsicherheit exp. uncertainty	Toleranz Tolerance	Konformität <sup>1)</sup> Conformity <sup>1)</sup>
1	500 g	500,01 g	0,01 g	0,016 g	0,050 g	✓
2	1000 g	1000,01 g	0,01 g	0,016 g	0,050 g	✓
3	1500 g	1500,01 g	0,01 g	0,016 g	0,050 g	✓
4	2000 g	2000,01 g	0,01 g	0,017 g	0,100 g	✓
5	3000 g	3000,02 g	0,02 g	0,018 g	0,100 g	✓

<sup>1)</sup> Bewertungskriterium: | [Abweichung] | + [erw. Unsicherheit] ≤ [Toleranz]  
Assessment criterion: | [Error] | + [exp. uncertainty] ≤ [Tolerance]

**Zusammenfassung / Summary**

Zum Zeitpunkt der Prüfung lagen die im Rahmen dieses Kalibrierscheins ermittelten Messergebnisse innerhalb der Toleranz.

At the time of testing, all measurement results determined in the context of this calibration certificate were within the tolerance.

Für einen anderen Verwendungsort, bei anderen Umgebungsbedingungen oder nicht bestimmungsgemäßer Verwendung der Waage kann Gültigkeit der oben aufgeführten Konformitätsaussagen nicht garantiert werden.

The validity of the statements of conformity cannot be guaranteed for different places of use, environmental conditions or improper use.

With a certificate of conformity you get a statement about whether the balance meets your defined requirements. In conjunction with a calibration certificate with accreditation symbol it serves as documented proof that the balance fulfils the required process demands. When doing this the process owner for the balance can select from different temperature specifications – depending on its individual requirements:

Conformity evaluation on the basis of the:	KERN		Price
Usage accuracy*	relative absolute	969-511 969-512	on request
Calibration results*	relative absolute	969-513 969-514	on request
Measurements as manufacturer or customer specification	other manuf. cust. spec. KERN device	969-515 969-516 969-517	on request on request 23,-

relative = % / absolute = g

\*as supplement to the calibration certificate with accreditation symbol  
(Details see [www.kern-lab.com](http://www.kern-lab.com))




# Calibration and verification prices for electronic balances

Accredited initial calibration and recalibration of balances at the KERN factory

Checking equipment	Initial calibration	Price € excl. of VAT ex works	Recalibration	Price € excl. of VAT ex works
Weighing range				
<b>Analytical balances</b>				
[Max] ≤ 5 kg	963-101	<b>210,-</b>	963-101R	<b>215,-</b>
[Max] > 5 kg	963-102	<b>260,-</b>	963-102R	<b>275,-</b>
<b>High resolution precision balances (&gt;500.000d)</b>				
[Max] ≤ 5 kg	963-103	<b>178,-</b>	963-103R	<b>183,-</b>
[Max] > 5 kg – 50 kg	963-104	<b>205,-</b>	963-104R	<b>210,-</b>
[Max] > 50 kg – 350 kg	963-105	<b>225,-</b>	963-105R	<b>235,-</b>
<b>Precision balances / industrial scales</b>				
[Max] ≤ 5 kg	963-127	<b>108,-</b>	963-127R	<b>111,-</b>
[Max] > 5 kg – 50 kg	963-128	<b>130,-</b>	963-128R	<b>133,-</b>
[Max] > 50 kg – 350 kg	963-129	<b>161,-</b>	963-129R	<b>165,-</b>
[Max] > 350 kg – 1500 kg	963-130	<b>225,-</b>	963-130R	<b>235,-</b>
[Max] > 1500 kg – 2900 kg <sup>1)</sup>	963-131	<b>305,-</b>	963-131R	<b>310,-</b>
[Max] > 2900 kg – 6000 kg <sup>1)</sup>	963-132	<b>610,-</b>	963-132R	<b>620,-</b>
[Max] > 6000 kg – 12000 kg <sup>1)</sup>	963-133	<b>680,-</b>	963-133R	<b>700,-</b>
<b>Hanging scales / crane scales</b>				
[Max] ≤ 5 kg	963-127H	<b>108,-</b>	963-127HR	<b>111,-</b>
[Max] > 5 kg – 50 kg	963-128H	<b>130,-</b>	963-128HR	<b>133,-</b>
[Max] > 50 kg – 350 kg	963-129H	<b>152,-</b>	963-129HR	<b>156,-</b>
[Max] > 350 kg – 1500 kg	963-130H	<b>270,-</b>	963-130HR	<b>285,-</b>
[Max] > 1500 kg – 2900 kg	963-131H	<b>415,-</b>	963-131HR	<b>420,-</b>
[Max] > 2900 kg – 6000 kg	963-132H	<b>680,-</b>	963-132HR	<b>700,-</b>
[Max] > 6000 kg – 12000 kg <sup>2) 3)</sup>	963-133H	<b>960,-</b>	963-133HR	<b>990,-</b>
<b>Additional services</b>				
Preparation for recalibration (cleaning, function test)			969-003R	<b>27,-</b>
Minimum weight of sample (for details see page 9 or internet)			969-103	<b>10,-</b>
Express service with delivery time 48 hours, per scale			962-116R	<b>55,-</b>
Additional test point during calibration			963-140R	<b>5,50</b>
Express shipping			in GER only	-

<sup>1)</sup> Floor scales & axle load scales only (Price per weighing panel). Please ask for further details. <sup>2)</sup> On request <sup>3)</sup> Processing time 4 working days

 Calibration prices for on-site calibration on request

## Verification<sup>6)</sup> prices for electronic balances

Checking equipment	Initial verification <sup>6)</sup>	Price € excl. of VAT ex works	Reverification <sup>7)</sup>	Price € excl. of VAT ex works
	<b>KERN</b>		<b>KERN</b>	
Electronic balances, class I, [Max] ≤ 5 kg <sup>4)</sup>	965-201	<b>158,-</b>	950-101R	<b>245,-</b>
Electronic balances, class I, [Max] > 5 kg <sup>4)</sup>	965-202	<b>158,-</b>	950-102R	<b>345,-</b>
Electronic balances, class II, [Max] ≤ 5 kg <sup>4)</sup>	965-216	<b>89,-</b>	950-116R	<b>136,-</b>
Electronic balances, class II, [Max] > 5 kg – 50 kg <sup>4)</sup>	965-217	<b>104,-</b>	950-117R	<b>164,-</b>
Electronic balances, class II, [Max] > 50 kg – 350 kg <sup>4)</sup>	965-218	<b>142,-</b>	950-118R	<b>235,-</b>
Electronic balances, class III-IV, [Max] ≤ 5 kg <sup>4)</sup>	965-227	<b>68,-</b>	950-127R	<b>120,-</b>
Electronic balances, class III-IV, [Max] > 5 kg – 50 kg <sup>4)</sup>	965-228	<b>88,-</b>	950-128R	<b>150,-</b>
Electronic balances, class III-IV, [Max] > 50 kg – 350 kg <sup>4)</sup>	965-229	<b>116,-</b>	950-129R	<b>205,-</b>
Electronic balances, class III-IV, [Max] > 350 kg – 1500 kg <sup>4)</sup>	965-230	<b>166,-</b>	950-130R	<b>345,-</b>
Electronic balances, class III-IV, [Max] > 1500 kg – 2900 kg <sup>4)</sup>	965-231	<b>187,-</b>	950-131R	<b>400,-</b>
Electronic balances, class III-IV, [Max] > 2900 kg – 6000 kg <sup>4)</sup>	965-232	<b>245,-</b>	950-132R	<b>610,-</b>
<b>Preparation for recalibration</b> (cleaning, function test)	-	-	969-006R	<b>26,-</b>

## Verification<sup>6)</sup> prices for electronic crane scales

Checking equipment	Initial verification <sup>6)</sup>	Price € excl. of VAT ex works	Reverification <sup>7)</sup>	Price € excl. of VAT ex works
	<b>KERN</b>		<b>KERN</b>	
Electronic crane scales, class III-IV, [Max] > 50 kg – 350 kg <sup>4)</sup>	965-129H	<b>125,-</b>	950-129HR	<b>210,-</b>
Electronic crane scales, class III-IV, [Max] > 350 kg – 1500 kg <sup>4)</sup>	965-130H	<b>158,-</b>	950-130HR	<b>460,-</b>
Electronic crane scales, class III-IV, [Max] > 1500 kg – 2900 kg <sup>4)</sup>	965-131H	<b>187,-</b>	950-131HR	<b>530,-</b>
Electronic crane scales, class III-IV, [Max] > 2900 kg – 6000 kg <sup>4)</sup>	965-132H	<b>245,-</b>	950-132HR	<b>760,-</b>
Electronic crane scales, class III-IV, [Max] > 6000 kg – 12000 kg <sup>4)</sup>	965-133H	<b>375,-</b>	950-133HR	<b>760,-</b>
Electronic crane scales, class III-IV, [Max] > 12000 kg – 31000 kg <sup>5)</sup>	-	-	950-134HR	<b>on request</b>
Electronic crane scales, class III-IV, [Max] > 31000 kg – 50000 kg <sup>5)</sup>	-	-	950-135HR	<b>on request</b>
<b>Preparation for recalibration</b> (cleaning, function test)	-	-	969-006R	<b>26,-</b>

<sup>4)</sup> Processing time 4 working days, <sup>5)</sup> processing time 15 working days,

<sup>6)</sup> "Initial verification"/conformity assessment according to NAWID: 2014/31/EU, only when purchasing a balance from KERN, valid throughout Europe

<sup>7)</sup> Verification ("reverification"), only for Germany

# Equipment qualification

Documented quality of your balances in the log book

Consistently high product quality requires the use of measuring and test equipment that provides comprehensible, consistent and reproducible results. Hence, quality management systems require that measuring and test equipment produces a detailed traceable description and documentation of calibration results and conformity statements. Work not documented is work not done.

Equipment qualification is documentary evidence that a equipment is suitable for the intended purpose and is working faultlessly. A balance log book as well as our EQS (Equipment qualification software) is used to record all activities and results required for the qualification and monitoring of balances during routine operation. This includes the installation and commissioning of the balances, routine tests, maintenance as well as the recording of special events (failures, repairs, change of location).

The structure of the balance log book is based on the qualification process of the balance. The requirements for the qualification system such as DIN EN ISO 9001, DIN EN ISO/IEC 17025, GLP/GMP, VDA must be taken into account. The log book supports the user in his/her daily work with the balance and is meant to serve as necessary evidence during inspections and audits. The responsibility for maintaining the log book and its appropriate use is to be borne by the user.

## Our proposal: Count on our support!

KERN offers this qualification concept throughout. Our validation services are carried out on the spot by technicians of our calibration laboratory and comprise among other things: installation, measurement test inclusive calibration certificate with accreditation symbol as well as records in your qualification log book and in the EQS (Equipment qualification software).

We give you advice about the options of device qualification, as required and will be happy to set up an appointment for qualification at the place of installation. We offer individual calibration and maintenance agreements for the periodically required requalification.

Further information can be found at [www.kern-lab.com](http://www.kern-lab.com)

**i** If you are interested in a qualification or training for equipment qualification, please feel free to contact us:  
**+49 7433 9933-400**  
**testservices-onsite@kern-sohn.com**

## Important elements of equipment qualification:



### Design Qualification (DQ)

With the design qualification, which is carried out under consideration of a requirement specification/functional description, all the requirements on which you as a user depend on, are defined. The purchase decision is made on the basis of the design specifications and the available devices. Careful selection in the DQ can prevent later deficiencies.



### Installation Qualification (IQ)

All steps to be taken for the installation and commissioning of the equipment are described in detail in the installation qualification. These include among others:

- checking for completeness of delivery and assurance that the delivered equipment meets the required specifications
- a description of the ambient conditions at the place of installation
- proper installation and assurance that the equipment is ready for operation after installation
- documentation of equipment configuration and equipment settings
- Recording and installation of connected peripherals units



### Function Qualification (OQ)

The operational qualification describes the metrological test performed for the balance at the place of installation. In the course of this all parameters that define the efficiency of a measurement will be checked. The OQ must be carried out by trained staff with the help of qualified aids (such as certified weights that are traceable to an approved standard). Briefing/training of users must be assured and recorded in the OQ.



### Performance Qualification (PQ)

The PQ represents documented evidence that the balance or weighing system functions in the selected application as intended. This will be assured by a qualification test of the equipment under real conditions with respect to its surroundings and the problem definition (such as traceable data transmission). If the balance or weighing system is "merely" to be used for weighing it will not be necessary to carry out a PQ as the ability to function has already been proven during the metrological test (OQ).



### Maintenance Qualification (MQ)

The periodical maintenance, cleaning work and complete metrological test of the balance/weighing system is documented in the MQ by a trained authorised engineer. The results are documented on a calibration certificate with accreditation symbol. Maintenance is carried out with the help of a maintenance schedule.

# The test weight

## Calibration of test weights

In order to have calibrated measuring devices you must have calibrated checking equipment. For balances, for example, the test equipment is calibrated test weights.

Depending on frequency of use, test weights must be recalibrated at regular intervals. This is the only way to guarantee that you maintain the requirements for checking equipment so that it functions reliably.

Recalibration times depend on the frequency of use, the conditions of use and your safety requirements.

The standard does not give a defined recalibration period. We recommend that, with intensive use, you recalibrate your test weights every 6 months and with normal use (daily), every 12 months (weekly).

### KERN calibrates test weights

- In all OIML error limit classes E1 – M3 and in sizes 1 mg – 2500 kg
- Test weights with free nominal value (any weight value)
- Carried out in Newton
- Independent of design (special designs)
- **Independent of the brand**



### The advantages of using KERN in-house calibration (you send your test weights to us)

- Excellent price to performance ratio
- Quick handling time
  - standard: approx. 4 working days
  - Express Service: from 48 hours (details on request)
- **Calibration service independent of the brand**
- KERN also reconditions old weights (e.g. cleaning or readjustment)
- KERN calibration certificates with accreditation symbol are valid internationally
- We would be pleased to monitor your recalibration times
- On request, collection and delivery service by our courier service
- The most modern calibration methods with robot operated comparators allow the most accurate calibration results and rapid process times



### The advantages of using KERN on-site calibration (we visit you)

We would be pleased to visit you within Germany and calibrate your test weights to OIML error limits M1 – M3, 10 kg – 2500 kg with our mobile MACOS calibration system. Minimized downtime of your checking equipment and direct contact with our expert are the major benefits of this service. Prices on request.



## Selection of the appropriate test weight

### Quality of the test weight

A balance can never be more accurate than the test weight which is used to adjust it. It depends on its tolerances.

### Accuracy of the test weight

Must approximately correspond to the readout [d] of the balance, or should if anything, be better.

### Weight size

This is often shown in the adjust mode "CAL" in the balance display. Given a choice, the largest displayed weight is the best one to use for accurate measurement. The weight of your test weight should ideally be larger than 80 % of the maximum weighing range of the balance. If accuracy and weight size (nominal value) are fixed, the appropriate test weight is selected according to the tolerances of the individual accuracy classes (error limit classes) E1 to M3 (see page 14).

### Example:

Balance with weighing range Max 2000 g (2 kg) and readout [d] 0.01 g (10 mg)

- The accuracy of the required test weight is determined by the readout [d] with approx.  $\pm 10$  mg.
- Displayed weight size in mode "CAL": 1000 g or 2000 g. The required test weight therefore has the weight size 2 kg.
- Appropriate test weight with tolerance  $\pm 10$  mg and weight size 2 kg is found in error limit class F1.

### Exception analytical balances (readout [d] $\leq 0.1$ mg):

E1 test weights are recommended. Depending on the safety requirement, E2 test weights with a calibration certificate with accreditation symbol will also be sufficient.

**For more information about our test weights please go to [kern-sohn.com](http://kern-sohn.com).**



# Standard OIML R 111 for weights

## The key points from the OIML norm R 111

OIML (Organisation Internationale de Metrologie Legale) has established the exact metrological requirements for weights in verified applications in approx. 100 states all over the world. The OIML recommendation R 111 for weights relates to sizes 1 mg – 5000 kg. Statements are made on the accuracy, materials, geometric shape, marking and storage of the weights.

## Error limits for weights of classes E1 to M3

The error limit classes are in fixed hierarchical levels in the proportion of 1:3, where E1 is the most accurate and M3 is the least accurate weight class. When testing weights with other weights, the correct test class is the next highest class.

## Error limit classes (= tolerances)

The values given in the table below (tolerances  $\pm \dots$  mg) are the respective permitted fabrication tolerances.

They are to be equal to the ► **measuring uncertainty** of the weight, if no ► **accredited calibration** is available.

## Conventional mass


The problem is the air buoyancy, which makes the weight appear lighter. In order to avoid this “distortion” in daily use, all weights are adjusted to the unit specifications as given in R 111, e.g. it is accepted that: material density of the weights is 8000 kg/m<sup>3</sup>, air density is 1.2 kg/m<sup>3</sup> and measuring temperature is 20 °C.

## KERN cylindrical test weights

Comply with OIML R 111 in all respects, without exception.

Nominal value ↓	OIML R 111:2004 Maximum permissible errors for weights = permissible tolerances “Tol ± mg”						
	E1	E2	F1	F2	M1	M2	M3
1 mg	± 0,003 mg	± 0,006 mg	± 0,020 mg	± 0,06 mg	± 0,20 mg	-	-
2 mg	± 0,003 mg	± 0,006 mg	± 0,020 mg	± 0,06 mg	± 0,20 mg	-	-
5 mg	± 0,003 mg	± 0,006 mg	± 0,020 mg	± 0,06 mg	± 0,20 mg	-	-
10 mg	± 0,003 mg	± 0,008 mg	± 0,025 mg	± 0,08 mg	± 0,25 mg	-	-
20 mg	± 0,003 mg	± 0,010 mg	± 0,03 mg	± 0,10 mg	± 0,3 mg	-	-
50 mg	± 0,004 mg	± 0,012 mg	± 0,04 mg	± 0,12 mg	± 0,4 mg	-	-
100 mg	± 0,005 mg	± 0,016 mg	± 0,05 mg	± 0,16 mg	± 0,5 mg	± 1,6 mg	-
200 mg	± 0,006 mg	± 0,020 mg	± 0,06 mg	± 0,20 mg	± 0,6 mg	± 2,0 mg	-
500 mg	± 0,008 mg	± 0,025 mg	± 0,08 mg	± 0,25 mg	± 0,8 mg	± 2,5 mg	-
1 g	± 0,010 mg	± 0,03 mg	± 0,10 mg	± 0,3 mg	± 1,0 mg	± 3,0 mg	± 10 mg
2 g	± 0,012 mg	± 0,04 mg	± 0,12 mg	± 0,4 mg	± 1,2 mg	± 4,0 mg	± 12 mg
5 g	± 0,016 mg	± 0,05 mg	± 0,16 mg	± 0,5 mg	± 1,6 mg	± 5,0 mg	± 16 mg
10 g	± 0,020 mg	± 0,06 mg	± 0,20 mg	± 0,6 mg	± 2,0 mg	± 6,0 mg	± 20 mg
20 g	± 0,025 mg	± 0,08 mg	± 0,25 mg	± 0,8 mg	± 2,5 mg	± 8,0 mg	± 25 mg
50 g	± 0,03 mg	± 0,10 mg	± 0,3 mg	± 1,0 mg	± 3,0 mg	± 10 mg	± 30 mg
100 g	± 0,05 mg	± 0,16 mg	± 0,5 mg	± 1,6 mg	± 5,0 mg	± 16 mg	± 50 mg
200 g	± 0,10 mg	± 0,3 mg	± 1,0 mg	± 3,0 mg	± 10 mg	± 30 mg	± 100 mg
500 g	± 0,25 mg	± 0,8 mg	± 2,5 mg	± 8,0 mg	± 25 mg	± 80 mg	± 250 mg
1 kg	± 0,5 mg	± 1,6 mg	± 5,0 mg	± 16 mg	± 50 mg	± 160 mg	± 500 mg
2 kg	± 1,0 mg	± 3,0 mg	± 10 mg	± 30 mg	± 100 mg	± 300 mg	± 1000 mg
5 kg	± 2,5 mg	± 8,0 mg	± 25 mg	± 80 mg	± 250 mg	± 800 mg	± 2500 mg
10 kg	± 5,0 mg	± 16 mg	± 50 mg	± 160 mg	± 500 mg	± 1600 mg	± 5000 mg
20 kg	± 10 mg	± 30 mg	± 100 mg	± 300 mg	± 1000 mg	± 3000 mg	± 10 g
50 kg	± 25 mg	± 80 mg	± 250 mg	± 800 mg	± 2500 mg	± 8000 mg	± 25 g
100 kg	-	± 160 mg	± 500 mg	± 1600 mg	± 5000 mg	± 16 g	± 50 g
200 kg	-	± 300 mg	± 1000 mg	± 3000 mg	± 10 g	± 30 g	± 100 g
500 kg	-	± 800 mg	± 2500 mg	± 8000 mg	± 25 g	± 80 g	± 250 g
1000 kg	-	± 1600 mg	± 5000 mg	± 16 g	± 50 g	± 160 g	± 500 g
2000 kg	-	-	± 10 g	± 30 g	± 100 g	± 300 g	± 1000 g
5000 kg	-	-	± 25 g	± 80 g	± 250 g	± 800 g	± 2500 g

## Composition table, valid for all KERN weight sets from 1 mg



# Calibration service for test weights (Accredited calibration)



The accreditation is valid for the scope defined in the document annex D-K-19408-01-00.

## Calibration certificate with accreditation symbol for weights

### 1. "Official" document

The calibration laboratory KERN (D-K-19408-01-00) is accredited through the accreditation point of the Deutsche Akkreditierungsstelle GmbH. The calibration certificate with accreditation symbol is recognised internationally and is available in several languages.

### 2. Item to be calibrated

The calibration item with nominal value and OIML tolerance class if applicable, as well as the serial number is documented. In this way the assignment of the issued calibration certificate to the weight or set of weights is completely guaranteed.

### 3. Traceability

The reference standards of the accredited laboratory are monitored in strictly defined cycles and periodically brought into line with national and thereby international standards. This is carefully documented and given on the calibration certificate. In this way the basic fundamental traceability to the national standard is ensured.

### 4. Applicant

On the very first page of the calibration certificate you will clearly see the applicant or owner of the calibrated checking equipment.

### 5. Environmental conditions

The environmental conditions during calibration are given here, such as the current temperature, the relative humidity and the air pressure at that point in time.

### 6. Metrological part

In this part of the calibration certificate information is given on the environmental conditions during calibration. Material, shape and density of the weight is given. The conventional weight value including the relevant measurement uncertainty is shown, as well as the OIML error limits and the OIML class.

### 7. ► Conventional mass

Using the substitution weighing method (comparative measurement with a test weight) you can determine the exact value of the weight to be calibrated. The conventional mass gives the deviation of the given value from the nominal value of the test item.

### 8. Measuring uncertainty

When obtaining any technical measurement, there is a particular uncertainty when trying to determine an exact value. This so-called measuring uncertainty should objectify measuring results, by establishing to what degree the measurement is expected to deviate from the true value. Determining and declaring the measuring uncertainty is of great significance, because the smaller this is, the more accurate the obtained weight.

**KERN** **KERN & SOHN GmbH**  
Akreditiertes Kalibrierlabor seit 1994.  
Accredited calibration laboratory since 1994.

Ihr Partner für Kalibrierdienstleistungen, Prüfmittelmanagement und Beratung.  
Your partner for calibration services, test equipment management and support.

Mitglied im / member of the  
**Deutschen Kalibrierdienst** **DKD** **ILAC-MRA** **DAKKS**  
Deutsche Akkreditierungsstelle  
D-K-19408-01-00

**Kalibrierschein** **Sample-2026-04/1** **Kalibrierzeichen** **Sample**  
Calibration certificate Calibration mark  
D-K-19408-01-00  
2026-04

Gegenstand **Gewichtssatz, 1 mg - 1 kg**  
Object **Set of weights, 1 mg - 1 kg**  
Klasse **E2**  
Class **E2**

Hersteller **KERN & Sohn GmbH**  
Manufacturer **Ziegele 1**  
**D-72336 Balingen**  
**Germany**

Typ **313-052**  
Type

Fabrikat/Serien-Nr. **G123456789**  
Serial number

Auftraggeber **Mustermann GmbH**  
Customer

Auftragsnummer **2026-123456789**  
Order No.

Anzahl der Seiten des Kalibrierscheines **3**  
Number of pages of the certificate

Kalibrierverfahren: **Die Kalibrierung erfolgte durch Vergleich mit den Bezugsnormalen des Kalibrierlaboratoriums nach der Substitutionsmethode mit Auftriebskorrektur.**  
Calibration method **The calibration ensued through comparison with the reference standards of the calibration laboratory using the substitution method with air buoyancy correction.**

Ort der Kalibrierung: **Kalibrierlaboratorium KERN**  
Place of calibration **Calibration - Laboratory KERN**

Umgebungsbedingungen: **Die Kalibrierung wurde bei folgenden Umgebungsbedingungen ausgeführt.**  
Ambient conditions **The calibration was carried out under the following ambient conditions:**

	von from	bis to	Unsicherheit uncertainty
Temperatur (°C) temperature	22,9	24,1	0,1
rel. Luftfeuchte (%) relative humidity	48,5	53,4	2,0
Luftdruck (hPa) air pressure	942,5	948,5	0,3

Magnetische Eigenschaften: **Der Hersteller hat bestätigt, dass die Gewichtsstücke die magnetischen Eigenschaften gemäß R111:2004 einhalten.**  
Magnetic properties **The manufacturer has confirmed the compliance of the magnetic properties of the weight pieces with the OIML R111:2004.**

Referenzgewichte: **123-D-K-19408-01-00-2022-05**  
Standard weights

Material / angenommene Dichte: **7850 kg/m³**  
Material / assumed density

Nennwert nominal value	Dichte density	Unsicherheit uncertainty
500 mg	7850 kg/m³	140 kg/m³

Messergebnisse: **2026-04**  
Measurement results

Nennwert nominal value	Kennzeichnung marking	konventioneller Wägewert conventional mass	Unsicherheit k=2 uncertainty	Fehlergrenze max. perm. error	Klasse <sup>a</sup> class <sup>a</sup>
1 mg		1 mg + 0.0010 mg	0.0020 mg	± 0.0080 mg	E2 ✓
2 mg		2 mg + 0.0005 mg	0.0020 mg	± 0.0080 mg	E2 ✓
2 mg	*	2 mg + 0.0016 mg	0.0020 mg	± 0.0080 mg	E2 ✓
5 mg		5 mg + 0.0010 mg	0.0020 mg	± 0.0080 mg	E2 ✓
10 mg		10 mg + 0.0009 mg	0.0020 mg	± 0.0080 mg	E2 ✓
20 mg		20 mg - 0.0001 mg	0.003 mg	± 0.010 mg	E2 ✓
20 mg	*	20 mg + 0.001 mg	0.003 mg	± 0.010 mg	E2 ✓
50 mg		50 mg + 0.001 mg	0.004 mg	± 0.012 mg	E2 ✓
100 mg		100 mg + 0.001 mg	0.005 mg	± 0.016 mg	E2 ✓
200 mg		200 mg + 0.002 mg	0.006 mg	± 0.020 mg	E2 ✓
200 mg	*	200 mg + 0.003 mg	0.006 mg	± 0.020 mg	E2 ✓
500 mg		500 mg + 0.005 mg	0.008 mg	± 0.025 mg	E2 ✓
1 g		1 g + 0.002 mg	0.010 mg	± 0.030 mg	E2 ✓
2 g		2 g + 0.002 mg	0.013 mg	± 0.040 mg	E2 ✓
2 g	*	2 g + 0.002 mg	0.013 mg	± 0.040 mg	E2 ✓
5 g		5 g + 0.010 mg	0.016 mg	± 0.050 mg	E2 ✓
10 g		10 g - 0.007 mg	0.020 mg	± 0.060 mg	E2 ✓
20 g		20 g + 0.005 mg	0.026 mg	± 0.080 mg	E2 ✓
20 g	*	20 g + 0.015 mg	0.026 mg	± 0.080 mg	E2 ✓
50 g		50 g + 0.02 mg	0.03 mg	± 0.10 mg	E2 ✓
100 g		100 g + 0.01 mg	0.05 mg	± 0.16 mg	E2 ✓
200 g		200 g + 0.05 mg	0.10 mg	± 0.30 mg	E2 ✓
200 g	*	200 g - 0.00 mg	0.10 mg	± 0.30 mg	E2 ✓
500 g		500 g + 0.10 mg	0.26 mg	± 0.80 mg	E2 ✓
1 kg		1 kg + 0.1 mg	0.5 mg	± 1.6 mg	E2 ✓

## Calibration certificate with accreditation symbol for test weights (extract)


Please see [www.kern-lab.com](http://www.kern-lab.com) for more details on our calibration service and other useful information

# Recalibration prices for test weights (Accredited calibration)

Class acc. OIML R 111:2004 →	E1 with volume determination (for new weights only)		E1 without volume determination		E2		F1 / F2 * F2 only		M1 / M2 / M3	
Nominal value ↓	KERN	Price € excl. of VAT ex works	KERN	Price € excl. of VAT ex works	KERN	Price € excl. of VAT ex works	KERN	Price € excl. of VAT ex works	KERN	Price € excl. of VAT ex works
1 mg	–	–	962-251R	87,-	962-351R	38,-	962-451R	24,-	962-651R	20,-
2 mg	–	–	962-252R	87,-	962-352R	38,-	962-452R	24,-	962-652R	20,-
5 mg	–	–	962-253R	87,-	962-353R	38,-	962-453R	24,-	962-653R	20,-
10 mg	–	–	962-254R	87,-	962-354R	38,-	962-454R	24,-	962-654R	20,-
20 mg	–	–	962-255R	87,-	962-355R	38,-	962-455R	24,-	962-655R	20,-
50 mg	–	–	962-256R	87,-	962-356R	38,-	962-456R	24,-	962-656R	20,-
100 mg	–	–	962-257R	87,-	962-357R	38,-	962-457R	24,-	962-657R	20,-
200 mg	–	–	962-258R	87,-	962-358R	38,-	962-458R	24,-	962-658R	20,-
500 mg	–	–	962-259R	87,-	962-359R	38,-	962-459R	24,-	962-659R	20,-
1 g	963-231	260,-	962-231R	87,-	962-331R	38,-	962-431R	24,-	962-631R	20,-
2 g	963-232	260,-	962-232R	87,-	962-332R	38,-	962-432R	24,-	962-632R	20,-
5 g	963-233	260,-	962-233R	87,-	962-333R	38,-	962-433R	24,-	962-633R	20,-
10 g	963-234	260,-	962-234R	87,-	962-334R	38,-	962-434R	24,-	962-634R	20,-
20 g	963-235	260,-	962-235R	87,-	962-335R	38,-	962-435R	24,-	962-635R	20,-
50 g	963-236	260,-	962-236R	87,-	962-336R	38,-	962-436R	24,-	962-636R	20,-
100 g	963-237	260,-	962-237R	87,-	962-337R	48,-	962-437R	27,-	962-637R	22,-
200 g	963-238	260,-	962-238R	87,-	962-338R	48,-	962-438R	27,-	962-638R	22,-
500 g	963-239	260,-	962-239R	87,-	962-339R	48,-	962-439R	27,-	962-639R	22,-
1 kg	963-241	260,-	962-241R	87,-	962-341R	48,-	962-441R	27,-	962-641R	22,-
2 kg	963-242	600,-	962-242R	107,-	962-342R	58,-	962-442R	35,-	962-642R	23,-
5 kg	963-243	600,-	962-243R	107,-	962-343R	58,-	962-443R	35,-	962-643R	23,-
10 kg	963-244	600,-	962-244R	107,-	962-344R	58,-	962-444R	35,-	962-644R	23,-
20 kg	963-245	1390,-	962-245R	810,-	962-345R	77,-	962-445R	40,-	962-645R	30,-
50 kg	963-246	1620,-	962-246R	900,-	962-346R	90,-	962-446R	55,-	962-646R	33,-
100 kg	–	–	–	–	–	–	962-591R*	160,-	962-691R	86,-
200 kg	–	–	–	–	–	–	962-592R*	160,-	962-692R	86,-
500 kg	–	–	–	–	–	–	962-593R*	160,-	962-693R	86,-
1000 kg	–	–	–	–	–	–	–	–	962-694R	188,-
2000 kg	–	–	–	–	–	–	–	–	962-695R	345,-
1 mg – 500 mg	–	–	962-250R	560,-	962-350R	260,-	962-450R	138,-	962-650R	82,-
1 mg – 50 g	963-201	1510,-	962-201R	920,-	962-301R	430,-	962-401R	225,-	962-601R	140,-
1 mg – 100 g	963-202	1650,-	962-202R	950,-	962-302R	470,-	962-402R	245,-	962-602R	146,-
1 mg – 200 g	963-203	1890,-	962-203R	1050,-	962-303R	540,-	962-403R	275,-	962-603R	164,-
1 mg – 500 g	963-204	2020,-	962-204R	1090,-	962-304R	580,-	962-404R	285,-	962-604R	172,-
1 mg – 1 kg	963-205	2140,-	962-205R	1180,-	962-305R	620,-	962-405R	300,-	962-605R	180,-
1 mg – 2 kg	963-206	2790,-	962-206R	1250,-	962-306R	680,-	962-406R	345,-	962-606R	198,-
1 mg – 5 kg	963-207	3130,-	962-207R	1300,-	962-307R	720,-	962-407R	360,-	962-607R	210,-
1 mg – 10 kg	963-208	3560,-	962-208R	1350,-	962-308R	780,-	962-408R	395,-	962-608R	215,-
1 g – 50 g	963-215	1090,-	962-215R	410,-	962-315R	177,-	962-415R	92,-	962-615R	54,-
1 g – 100 g	963-216	1190,-	962-216R	445,-	962-316R	210,-	962-416R	106,-	962-616R	65,-
1 g – 200 g	963-217	1450,-	962-217R	530,-	962-317R	280,-	962-417R	134,-	962-617R	80,-
1 g – 500 g	963-218	1580,-	962-218R	580,-	962-318R	320,-	962-418R	150,-	962-618R	90,-
1 g – 1 kg	963-219	1730,-	962-219R	630,-	962-319R	355,-	962-419R	165,-	962-619R	96,-
1 g – 2 kg	963-220	2430,-	962-220R	720,-	962-320R	440,-	962-420R	205,-	962-620R	117,-
1 g – 5 kg	963-221	2840,-	962-221R	740,-	962-321R	495,-	962-421R	225,-	962-621R	126,-
1 g – 10 kg	963-222	3310,-	962-222R	810,-	962-322R	540,-	962-422R	250,-	962-622R	136,-

Additional costs for preparation, overhaul and adjustment before the calibration	KERN	Price € excl. of VAT ex works
<b>Preparation of weights (e.g. cleaning, etc.)</b>		
Single weight	969-001R	5,-
Weight set	969-002R	22,-
<b>Subsequent services are carried out after confirmation</b>		
Continued overhaul of weights (e.g. wet-cleaning, markings, repair, special packaging, adjustment E1, E2 ...)	969-005R	<b>T &amp; M basis</b>
Adjustment, per weight only available for weights with adjustment chamber (F1–M3)	969-010R	16,-
<b>Second calibration after adjustment or substitution, per weight</b>		
Class E1	969-210R	67,-
Class E1 incl. volume determination	969-211R	245,-
Class E2	969-310R	32,-
Class F1 / F2	969-410R	21,-
Class M1–M3	969-610R	17,-
<b>Testing of magnetic properties according to OIML R 111-2004, per weight</b>	961-115R	17,-
<b>Calibration of NON-OIML test weights, additional price per weight</b>	–	8,-

KERN delivery times	
<b>Standard service</b> Class E2–M3	4 working days
<b>Standard service</b> Class E1, 1 mg – 500 mg, and recalibration 1 g – 10 kg with a known volume	10 working days
Class E1, 1 g – 2 kg, incl. volume determination (new weights)	15 working days

	<b>Express service in 48 hours</b> except for class E1
<ul style="list-style-type: none"> <li>Urgent order is received at KERN by 12:00 noon at the latest</li> <li>Ready for shipping at KERN within two working days, at 12:00 noon</li> <li>Return by standard parcel service or express shipping (Costs and processing time on request)</li> <li>Additional cost for Express Service, for each</li> <li>KERN test weight KERN 962-115R € 23,-</li> <li>For Express shipping (details on request)</li> </ul>	



## Reverification prices for test weights

Class acc. OIML R 111 →	E2 with verification certificate		F1 / F2 with verification certificate		M1 with verification certificate	
Nominal value ↓	KERN	Price € excl. of VAT ex works	KERN	Price € excl. of VAT ex works	KERN	Price € excl. of VAT ex works
1 mg	952-351R	66,-	952-451R	57,-	952-651R	38,-
2 mg	952-352R	66,-	952-452R	57,-	952-652R	38,-
5 mg	952-353R	66,-	952-453R	57,-	952-653R	38,-
10 mg	952-354R	66,-	952-454R	57,-	952-654R	38,-
20 mg	952-355R	66,-	952-455R	57,-	952-655R	38,-
50 mg	952-356R	66,-	952-456R	57,-	952-656R	38,-
100 mg	952-357R	66,-	952-457R	57,-	952-657R	38,-
200 mg	952-358R	66,-	952-458R	57,-	952-658R	38,-
500 mg	952-359R	66,-	952-459R	57,-	952-659R	38,-
1 g	952-331R	66,-	952-431R	57,-	952-631R	38,-
2 g	952-332R	66,-	952-432R	57,-	952-632R	38,-
5 g	952-333R	66,-	952-433R	57,-	952-633R	38,-
10 g	952-334R	66,-	952-434R	57,-	952-634R	38,-
20 g	952-335R	66,-	952-435R	57,-	952-635R	38,-
50 g	952-336R	66,-	952-436R	57,-	952-636R	38,-
100 g	952-337R	73,-	952-437R	59,-	952-637R	39,-
200 g	952-338R	73,-	952-438R	59,-	952-638R	39,-
500 g	952-339R	73,-	952-439R	59,-	952-639R	39,-
1 kg	952-341R	73,-	952-441R	59,-	952-641R	39,-
2 kg	952-342R	91,-	952-442R	64,-	952-642R	40,-
5 kg	952-343R	91,-	952-443R	64,-	952-643R	40,-
10 kg	952-344R	91,-	952-444R	64,-	952-644R	52,-
20 kg	952-345R	92,-	952-445R	67,-	952-645R	57,-
50 kg	952-346R		952-446R	77,-	952-646R	58,-
1 mg – 500 mg	952-350R	435,-	952-450R	225,-	952-650R	141,-
1 mg – 50 g	952-301R	700,-	952-401R	370,-	952-601R	235,-
1 mg – 100 g	952-302R	760,-	952-402R	405,-	952-602R	250,-
1 mg – 200 g	952-303R	870,-	952-403R	450,-	952-603R	280,-
1 mg – 500 g	952-304R	930,-	952-404R	465,-	952-604R	290,-
1 mg – 1 kg	952-305R	960,-	952-405R	490,-	952-605R	310,-
1 mg – 2 kg	952-306R	1120,-	952-406R	560,-	952-606R	340,-
1 mg – 5 kg	952-307R	1200,-	952-407R	600,-	952-607R	360,-
1 mg – 10 kg	952-308R	1260,-	952-408R	640,-	952-608R	370,-
1 g – 50 g	952-315R	285,-	952-415R	164,-	952-615R	106,-
1 g – 100 g	952-316R	340,-	952-416R	174,-	952-616R	116,-
1 g – 200 g	952-317R	440,-	952-417R	220,-	952-617R	136,-
1 g – 500 g	952-318R	500,-	952-418R	245,-	952-618R	152,-
1 g – 1 kg	952-319R	560,-	952-419R	265,-	952-619R	166,-
1 g – 2 kg	952-320R	680,-	952-420R	340,-	952-620R	198,-
1 g – 5 kg	952-321R	750,-	952-421R	370,-	952-621R	220,-
1 g – 10 kg	952-322R	840,-	952-422R	410,-	952-622R	230,-

! Reverification only valid in Germany

KERN reverification delivery time	
Standard reverification service Class E2 – M1	6 working days

Additional costs for preparation, overhaul and adjustment before the reverification	KERN	Price € excl. of VAT ex works
---	------	--

Preparation of weights (e.g. cleaning, etc.)		
Single weight	969-008R	5,-
Weight set	969-009R	22,-

### Subsequent services are carried out after confirmation

Continued overhaul of weights (e.g. wet-cleaning, markings, repair, special packaging, adjustment E2)	969-005R	at cost
--	----------	---------

Adjustment, per weight only available for weights with adjustment chamber (F-M1)	969-010R	16,-
--	----------	------

### Verification after adjustment or substitution, per weight

Class E2	969-310R	32,-
Class F1 / F2	969-410R	21,-
Class M1	969-610R	17,-



# The force gauge

Accredited calibration with calibration certificate for force gauges

The KERN calibration laboratory is at your side when you need to calibrate accredited reliably.

From the transducer to the full measuring chain, we are happy to take care of traceable calibration of your test equipment for you. Our accreditation includes the calibration of tensile and pressure force up to 5 kN according to the standards DIN EN ISO 376 and DKD-R 3-3, each with the Newton (N) display unit for a complete measuring chain (situation A) or voltage ratio/transmission coefficient (mV/V, situation B).

Below you will find a comparison of which standard meets which criteria:

Comparison of DIN EN ISO 376 and DKD-R 3-3		
	ISO 376	DKD-R 3-3
Standard	ISO standard (internationally standardized)	Standard of the DKD (Germany)
Measuring equipment	Force transducers and complete measuring chains	Force transducers and complete measuring chains
Area of application	Specifically force gauges for the testing of testing equipment	General force gauges
Number of power stages	8	5
Classification/Assessment	Classification in classes 00; 0,5; 1 and 2	None in standard
Test sequences	Fixed procedure	Sequences A, B, C, D possible Standard is sequence A B, C and D are reduced sequences, relevant previous knowledge is necessary
Summary	Higher-quality calibration, as 8 force levels are calibrated	High-quality calibration, reduced sequences with less effort possible

We can offer you a calibration solution for the following situations:



**Situation A:**  
separate force transducer,  
display unit **mV/V**

**Situation B:**  
complete force gauge  
(consisting of transducer,  
amplifier and display),  
display unit **N**

You can find further information on this topic at:

**www.kern-lab.com**

**KERN & SOHN GmbH**  
Akkreditiertes Kalibrierlabor seit 1994.  
Accredited calibration laboratory since 1994.  
Ihr Partner für Kalibrierdienstleistungen, Prüfmittelmanagement und Beratung.  
Your partner for calibration services, test equipment management and support.

Mitglied im / member of the  
**Deutschen Kalibrierdienst** **DKD**   
DAkkS ist Unterzeichner der multilateralen Übereinkommen der European co-operation for Accreditation (EA) und der International Laboratory Accreditation Cooperation (ILAC) zur gegenseitigen Anerkennung der Kalibrierzertifikate. Für die Einhaltung einer angemessenen Frist zur Wiederholung der Kalibrierung ist der Benutzer verantwortlich.

**Kalibrierschein**  
Calibration Certificate

**Sample-2026-01/1**

**Kalibrierzeichen**  
Calibration mark

Sample

D-K-19408-01-00

2026-01

Gegenstand Object	Kraftmessgerät Force gauge  Max. 1000 N, d = 0,5 N	Dieser Kalibrierschein dokumentiert die metrologische Rückführbarkeit auf nationale Normale zur Darstellung der Einheiten in Übereinstimmung mit dem internationalen Einheitensystem (SI). Die DAkkS ist Unterzeichner der multilateralen Übereinkommen der European co-operation for Accreditation (EA) und der International Laboratory Accreditation Cooperation (ILAC) zur gegenseitigen Anerkennung der Kalibrierzertifikate. Für die Einhaltung einer angemessenen Frist zur Wiederholung der Kalibrierung ist der Benutzer verantwortlich.
Hersteller Manufacturer	Sauter GmbH Ziegel 1 72336 Balingen Deutschland	
Typ Type	FH 1K	
Seriennummer Serial number	SA20H02287	
Auftraggeber Customer	Musterfirma GmbH Musterstraße 1	

**Messwerte (Zugkraft) / Measurement results (tension force)**

Ausrichtung rotation	Ausgangsposition / initial position 0°		120°		240°	
Kraft force	R1	R2	R3	R4	R5	R6
0 N	0,0 N	0,0 N	0,0 N	0,0 N	0,0 N	0,0 N
200,0 N	-199,5 N	-199,5 N	-199,5 N	-200,0 N	-199,5 N	-200,0 N
400,0 N	-399,5 N	-399,5 N	-399,5 N	-399,5 N	-399,5 N	-399,5 N
600,0 N	-599,5 N	-599,5 N	-599,5 N	-599,5 N	-599,5 N	-599,5 N
800,0 N	-799,0 N	-799,0 N	-799,0 N	-799,0 N	-799,0 N	-799,0 N
1000,0 N	-999,5 N	-999,0 N	-999,0 N	-999,0 N	-999,0 N	-999,0 N
0 N	0,0 N	0,0 N	0,0 N	0,0 N	0,0 N	0,0 N

**Messergebnisse (Zugkraft) / Measured values (tension force)**  
Aus den oben aufgeführten Messwerten ergeben sich die folgenden Messergebnisse:  
The following measurement results are calculated using the measured values above:  
**Rel. Kalibrierwertabweichung / Rel. cal. max. load error:**  $b_1 = 0,000\%$   
**Rel. Nullpunktabweichungen / Rel. zero error:**  $f_0 = 0,000\%$  (R1), 0,000 % (R2), 0,000 % (R3/R4), 0,000 % (R5/R6)

Kraft force	arith. Mittelwert $\bar{x}$ average $\bar{x}$	rel. Wiederholpräzision $\bar{v}$ rel. repeatability $\bar{v}$	rel. Vergleichspräzision $\bar{v}$ rel. reproducibility $\bar{v}$	rel. Unsicherheitsgrenze $\bar{v}$ system $\bar{v}$
200,0 N	-199,5 N	0,000 %	0,000 %	+0,251 %
400,0 N	-399,5 N	0,000 %	0,000 %	0,000 %
600,0 N	-599,5 N	0,000 %	0,000 %	0,000 %
800,0 N	-799,0 N	0,000 %	0,000 %	+0,003 %
1000,0 N	-999,5 N	+0,000 %	+0,000 %	0,000 %

Calibration certificate with accreditation symbol for  
force gauges (extract)

## Prices for accredited recalibration of force gauges and force transducers

Situation A: Force transducer (voltage ratio, in mV/V)*1,2					
ISO 376 (8 stages)			DKD-R 3-3 (5 stages, sequence A)		
KERN	Measuring range	€	KERN	Measuring range	€
Tensile force:					
963-161IVR	≤ 500 N	270,-	963-161VR	≤ 500 N	255,-
963-162IVR	≤ 2 kN	325,-	963-162VR	≤ 2 kN	300,-
963-163IVR	≤ 5 kN	420,-	963-163VR	≤ 5 kN	390,-
Compression force:					
963-261IVR	≤ 500 N	270,-	963-261VR	≤ 500 N	255,-
963-262IVR	≤ 2 kN	325,-	963-262VR	≤ 2 kN	300,-
963-263IVR	≤ 5 kN	420,-	963-263VR	≤ 5 kN	390,-
Tensile & Compression force:					
963-361IVR	≤ 500 N	455,-	963-361VR	≤ 500 N	420,-
963-362IVR	≤ 2 kN	540,-	963-362VR	≤ 2 kN	500,-
963-363IVR	≤ 5 kN	720,-	963-363VR	≤ 5 kN	660,-
969-003R	Preparation for recalibration (cleaning, function test)				27,-

Situation B: Complete force gauge (in N)*2					
ISO 376 (8 stages)			DKD-R 3-3 (5 stages, sequence A)		
KERN	Measuring range	€	KERN	Measuring range	€
Tensile force:					
963-161IR	≤ 500 N	220,-	963-161R	≤ 500 N	200,-
963-162IR	≤ 2 kN	270,-	963-162R	≤ 2 kN	245,-
963-163IR	≤ 5 kN	375,-	963-163R	≤ 5 kN	340,-
Compression force:					
963-261IR	≤ 500 N	220,-	963-261R	≤ 500 N	200,-
963-262IR	≤ 2 kN	270,-	963-262R	≤ 2 kN	245,-
963-263IR	≤ 5 kN	375,-	963-263R	≤ 5 kN	340,-
Tensile & Compression force:					
963-361IR	≤ 500 N	305,-	963-361R	≤ 500 N	365,-
963-362IR	≤ 2 kN	495,-	963-362R	≤ 2 kN	455,-
963-363IR	≤ 5 kN	670,-	963-363R	≤ 5 kN	600,-
969-003R	Preparation for recalibration (cleaning, function test)				27,-

For each force gauge without interface or from other manufacturers we charge a surcharge for the additional effort.

\*1 Compatibility with our amplifiers required

\*2 Installation in our measuring equipment required



# Factory calibration

for force

Situation A: Force transducer (voltage ratio, in mV/V)*1,2			Situation B: Complete force gauge (in N)*2		
KERN	Measuring range	€	KERN	Measuring range	€
<b>Tensile force:</b>					
961-161VR	≤ 500 N	<b>255,-</b>	961-161R	≤ 500 N	<b>200,-</b>
961-162VR	≤ 2 kN	<b>300,-</b>	961-162R	≤ 2 kN	<b>245,-</b>
961-163VR	≤ 5 kN	<b>390,-</b>	961-163R	≤ 5 kN	<b>340,-</b>
961-164VR	≤ 20 kN	<b>495,-</b>	961-164R	≤ 20 kN	<b>445,-</b>
961-165VR	≤ 50 kN	<b>495,-</b>	961-165R	≤ 50 kN	<b>445,-</b>
961-166VR	≤ 120 kN	<b>530,-</b>	961-166R	≤ 120 kN	<b>490,-</b>
961-167VR	≤ 250 kN	<b>530,-</b>	961-167R	≤ 250 kN	<b>490,-</b>
<b>Compression force:</b>					
961-261VR	≤ 500 N	<b>255,-</b>	961-261R	≤ 500 N	<b>200,-</b>
961-262VR	≤ 2 kN	<b>300,-</b>	961-262R	≤ 2 kN	<b>245,-</b>
961-263VR	≤ 5 kN	<b>390,-</b>	961-263R	≤ 5 kN	<b>340,-</b>
961-264VR	≤ 20 kN	<b>495,-</b>	961-264R	≤ 20 kN	<b>445,-</b>
961-265VR	≤ 50 kN	<b>495,-</b>	961-265R	≤ 50 kN	<b>445,-</b>
961-266VR	≤ 120 kN	<b>530,-</b>	961-266R	≤ 120 kN	<b>490,-</b>
961-267VR	≤ 250 kN	<b>530,-</b>	961-267R	≤ 250 kN	<b>490,-</b>
<b>Tensile &amp; Compression force</b>					
961-361VR	≤ 500 N	<b>420,-</b>	961-361R	≤ 500 N	<b>365,-</b>
961-362VR	≤ 2 kN	<b>500,-</b>	961-362R	≤ 2 kN	<b>455,-</b>
961-363VR	≤ 5 kN	<b>660,-</b>	961-363R	≤ 5 kN	<b>600,-</b>
961-364VR	≤ 20 kN	<b>710,-</b>	961-364R	≤ 20 kN	<b>660,-</b>
961-365VR	≤ 50 kN	<b>710,-</b>	961-365R	≤ 50 kN	<b>660,-</b>
961-366VR	≤ 120 kN	<b>780,-</b>	961-366R	≤ 120 kN	<b>720,-</b>
961-367VR	≤ 250 kN	<b>780,-</b>	961-366R	≤ 250 kN	<b>720,-</b>
969-003R	Preparation for recalibration (cleaning, function test)				<b>27,-</b>

For each force gauge without interface or from other manufacturers we charge a surcharge for the additional effort.

\*1 Compatibility with our amplifiers required

\*2 Installation in our measuring equipment required

## The temperature and relative humidity

Accredited calibration with calibration certificate for temperature and relative humidity

We perform accredited calibrations according to DKD-R 5-1 and DKD-R 5-8 for measuring instruments for the recording of ambient conditions. Our scope of accreditation covers a measuring range of 5°C to 50°C for temperature sensors and a measuring range from 20 % to 75 % relative humidity for humidity sensors.


Checking equipment	KERN	Price € excl. of VAT ex works
Temperature measuring device, external sensor	963-613R	<b>143,-</b>
Temperature measuring device, internal sensor	963-623R	<b>143,-</b>
Temperature and humidity, combi-sensor, external sensor, 1 temperature & 3 humidity points	963-631R	<b>220,-</b>
Temperature and humidity, combi-sensor, external sensor, 3 temperature & 3 humidity points	963-633R	<b>360,-</b>
Temperature and humidity, combi-sensor, internal sensor, 1 temperature & 3 humidity points	963-641R	<b>220,-</b>
Temperature and humidity, combi-sensor, internal sensor, 3 temperature & 3 humidity points	963-643R	<b>360,-</b>
Temperature - additional test point	963-605R	<b>25,-</b>
Humidity - additional test point	963-606R	<b>25,-</b>
<b>Additional services</b>		
Preparation for recalibration (cleaning, functional test)	969-003R	<b>27,-</b>

For each measuring device without interface we charge a surcharge for the additional effort.



# Factory calibration

for other measuring instruments



**Factory calibration certificate for wrench testing devices**  
(extract). Further details on the internet [www.kern-lab.com](http://www.kern-lab.com)

## Factory calibration certificates

As calibration certificates with accreditation symbol cannot be offered for all measuring devices or measurement sizes, or where it is not customary, we then offer factory calibration certificates. These calibration certificates meet international standards and are particularly suitable as proof of exacting calibration in the monitoring of your checking equipment. Factory calibrations are available for many measuring devices, for example

- Mechanical balances (spring balances, etc.)
- Force-measuring devices up to 120 kN
- Measuring devices for layer thickness 0 µm – 2000 µm
- Hardness testing devices in accordance with Leeb tests
- Ultrasound material thickness testing devices 25 – 300 mm

This is not an accredited calibration (no proof of metrological traceability).

**We carry out calibrations whatever the brand.** In order to avoid any unnecessary delays when processing your order, please send us the technical documents and accessories with the checking device. Calibration time 4 working days.

KERN	Measurand	Measuring range	Price € excl. of VAT ex works
Factory calibration			
961-102KR	Force (for hand grip dynamometer KERN MAP)	≤ 130 kg	<b>180,-</b>
961-110R	Coating thickness	≤ 2000 µm F or N	<b>180,-</b>
961-112R	Coating thickness	≤ 2000 µm FN	<b>255,-</b>
961-113R	Wall thickness (ultra sound)	≤ 300 mm (in stainless steel)	<b>180,-</b>
961-170R	Hardness Shore	For sets up to 7 plates	<b>143,-</b>
961-131R	Hardness Leeb	400 – 800 HLD	<b>180,-</b>
961-132R	Hardness Leeb	Test block (for Leeb durometer)	<b>180,-</b>
961-270R	Hardness (UCI)	200 – 800 HV	<b>390,-</b>
961-150R	Length	≤ 300 mm	<b>180,-</b>
961-190R	Light	≤ 200000 lx	<b>350,-</b>
961-100R	Weight (Mechanical balances/ Spring balances)	≤ 5 kg	<b>107,-</b>
961-101R	Weight (Mechanical balances/ Spring balances)	> 5 – 50 kg	<b>133,-</b>
961-102R	Weight (Mechanical balances/ Spring balances)	> 50 – 350 kg	<b>158,-</b>
961-103R	Weight (Mechanical balances/ Spring balances)	> 350 – 1500 kg	<b>245,-</b>
961-104R	Weight (Mechanical balances/ Spring balances)	> 1500 – 3000 kg	<b>330,-</b>
961-105R	Weight (Mechanical balances/ Spring balances)	> 3000 – 6000 kg	<b>660,-</b>
961-106R	Weight (Mechanical balances/ Spring balances)	> 6000 – 12000 kg	<b>750,-</b>
961-120R	Wrench testing devices	1 Nm – 200 Nm	<b>255,-</b>
961-290R	Refractometer		<b>146,-</b>
964-305R	Temperature calibration moisture analyser		<b>183,-</b>
Additional services			
962-116R	Calibration express service with 48 hour delivery (only on new purchases)		<b>55,-/</b> Instrument
969-003R	Preparation for recalibration (cleaning, function test)		<b>27,-</b>

\* See [www.kern-lab.com](http://www.kern-lab.com) for an overview of calibratable instruments and test services for further measuring sizes.

# Digital calibration certificate (DCC)

The calibration certificate with accreditation symbol becomes digital

The digital age is constantly revealing new opportunities for technological innovations. As a result, topics such as sustainability and resource-saving work are becoming increasingly important, especially in the industrial context. Calibration certificates, as metrological proof of a check of the respective measurand, are still printed on paper and sent to the customer by post. However, the option of a calibration certificate in PDF format is already available. Both variants no longer meet the current requirements of a laboratory. In cooperation of the PTB (Physikalisch Technische Bundesanstalt) and several expert groups from industry, the idea of a digital calibration certificate (DCC) was born. Based on the file format XML, this should be both human- and machine-readable. The basic idea is to transmit the information contained in a conventional calibration certificate to the customer electronically and securely as data. KERN DCCs comply with the requirements of VDI/VDE 2623 and PTB.

- You can download your digital calibration certificate (DCC) on [www.kern-lab.com/dcc](http://www.kern-lab.com/dcc)

- This service is **free of charge!**

- The Digital Calibration Certificate (DCC) is currently available for the calibration of weights. Gradually, this will also be available for other measuring devices.



For classical calibration certificates there is a scheme predefined by standards, to which the calibration laboratories shall adhere in order to be or remain accredited. This ensures standardization and simplification of the calibration procedure and its documentation. PTB is constantly working on designing such a scheme also for the DCC. Within such a scheme, of course, all previous information of a calibration certificate shall be included. This is implemented by using a so-called XML schema file in XSD format. In conclusion the generated XML file can be checked against the schema, which guarantees that the DCC is well-formed and complies with the specifications of the PTB.

Appropriate encryption algorithms and a digital signature ensure that no one can subsequently change the document and customers have the security of knowing that their tested measuring equipment has been processed and documented according to maximum quality requirements.

Within this context, KERN not only stands for precision and quality, but also for innovation in all technical matters.

## YOUR ADVANTAGES OF THE DIGITAL CALIBRATION CERTIFICATE:

- ✓ Meets the requirements of DIN EN ISO/IEC 17025:2018
- ✓ Secure and error-free data transmission
- ✓ Automatic transfer of data into digital processes
- ✓ Evaluation of measurement data and optimisation of measurements
- ✓ Forgery protection through digital signature
- ✓ Internationally recognised through XML-based format
- ✓ Machine readability, supporting the quality monitoring processes

# KERN GLOSSARY

## CALIBRATION

Calibration is the testing and determination of the precision of a measure value without intervention in the measurement system. The calibration certificate contains the measured value with information on the relevant measuring uncertainty. If applicable, a statement can be made as to whether this is within tolerance limits. Industry requires calibration of measuring devices, in order to, for example, be able to connect parts manufactured at different locations without encountering problems. Calibrations must be repeated at appropriate time intervals, for which the user is responsible. KERN recommends that, with intensive (daily) use, you recalibrate the measuring devices every 6 months and with normal (weekly) use, every 12 months.

## ACCREDITED CALIBRATION

Accredited calibration is carried out for measuring devices, reference materials and material measures for particular measurement sizes and measurement ranges, which are defined individually for every laboratory as part of their accreditation. The issued calibration certificates with accreditation symbol are proof of the metrological traceability to national and international standards, as required, for example, by the DIN EN ISO 9000 and DIN EN ISO/IEC 17025 standards. Accredited calibration has no legally regulated period of validity. The operator is responsible for observing an appropriate time for recalibration.

## INTERNATIONAL VALIDITY OF CALIBRATION CERTIFICATES WITH ACCREDITATION SYMBOL

DAkkS is represented in the EA (European co-operation for Accreditation) as well as in the ILAC (International Laboratory Accreditation Cooperation). This ensures that accredited calibration and calibration certificates with accreditation symbol are recognised and valid almost anywhere in the world.

## ADJUSTING

Precise setting of a measuring device by professional intervention in the measurement system. For balances: Either with an external test weight using the adjustment function (CAL or CAL key), or with the automatic internal adjustment or adjustment control. This is necessary following changes in temperature, changed environmental conditions, change of location, etc. daily routine checks are recommended. The term "calibrating" was formerly also used for adjusting, but today it means something else (see above).

## MONITORING YOUR CHECKING EQUIPMENT

This is a mandatory requirement of quality management systems.

## TRACEABILITY

The precondition of a perfect measurement is the complete proof that a measuring device is traceable to the SI units. This is also a requirement of the most important international standards. For example the correct display of balances and test weights is checked and – if necessary – corrected (adjusted). The test weights are traced back through a chain of calibrations to the national standard at the PTB (Physikalisch Technische Bundesanstalt) in Braunschweig, which in turn is traceable to the SI definition via various physical realizations (the "watt balance" or the "silicon sphere").

The correct traceability is what makes international measurements possible in the first place.

## MEASURING UNCERTAINTY

Measurement uncertainty is determined for each balance according to a precisely given test method and documented in the Calibration certificate. It depends on various factors, both internal and external. The measuring uncertainty of a measuring device is an objective measure of its accuracy and is therefore an accurate statement for its appropriate use.

## OIML

OIML (Organisation Internationale de Métrologie Légale) has representatives from almost 100 states who work on unified build and test regulations for all measuring devices. In the OIML certification system the certificates issued by the member states certify that a particular measuring device build type is in accordance with the OIML recommendations. In this way, a build type which was tested and approved in one country, can be approved in another country without having to repeat the test. (Excerpt from PTB). The OIML R111 guideline defines the construction-related characteristics for test weights, such as material, upper surface texture, markings, construction, shape etc.

## FACTORY CALIBRATION CERTIFICATES

The testing of measuring devices for accuracy in accordance with a recognised but not accredited process without proof of metrological traceability – this is the difference when compared with accredited calibration.


## CONVENTIONAL MASS

Every body experiences a relatively small loss of weight in air (buoyancy). This must be taken into account for accurate weighing procedures. In order to avoid this "distortion" in daily use, all weights are adjusted to the unit specifications as given in R111 OIML recommendation. (air pressure  $1.2 \text{ kg/m}^3$  and material density  $8000 \text{ kg/m}^3$ )

**KERN & SOHN GmbH**  
 Calibration Service  
 Ziegelei 1  
 72336 Balingen  
 Germany  
 Tel. +49 7433 9933-196  
 Fax +49 7433 9933-149  
 testservices@kern-sohn.com  
 www.kern-sohn.com  
 www.kern-lab.com

## Accredited Calibration Laboratory for more than 30 years

QM certification and accreditation by KERN  
as a basis for the highest level of quality.



DIN EN ISO/IEC 17025:2018



NAWI: 2014/31/EU

**www.kern-lab.com –**

## The Central Portal for everything you need to know about the extensive KERN Calibration Services

On our website you will always find the latest news and useful information about testing and measuring devices, calibration, legal metrology and expansions to our range of services. You will also find numerous online services on the website

### Database supported management of test equipment

Information on your test equipment which has been calibrated by us is stored in our database. In this way it is possible to make trend calculations. You will therefore get an overview of the long-term stability and trend behaviour of your test equipment and the necessary recalibration period can easily be determined and specified.

### Paperless documentation

So that there is no administrative effort, we can handle all calibration documentation in a paperless process. From quotation, through to order confirmation, delivery note and invoice right up to calibration certificate, you will receive all documents by e-mail or you can retrieve them online.

Would you prefer to receive your calibration certificate or your invoice, for example, in paper form? Of course this is not a problem either.

### Price quote generator

Create your own offer – you will receive your offer directly and without delay.

### RMA (Return Material Authorization)

Using the quotation generator, you can have a Return Material Authorization (RMA) number created directly for sending your test equipment. This makes it very easy to send in your test equipment and to start the calibration directly after arrival in our laboratory!

This brochure is valid until a new version of the brochure is released. In Europe, all prices do not include the applicable V.A.T. You can find our terms & conditions at [www.kern-lab.com](http://www.kern-lab.com). The logo used is a registered trademark of Deutsche Akkreditierungsstelle (DAkkS).

